REF 505-11 **IVD** A usage diagnostique in vitro

Usage prévu et application

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Bibliographie

- Calatzis A, Spannagl M, Gempeler-Messina P, Kolde HJ, Schramm W, Haas S. The prothrombinase induced clotting test: A new technique for the monitoring of anticoagulants. Haemostasis 2000; 30 (Suppl. 2): 172-174
- National Committee for Clinical Laboratory Standards. Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays. NCCLS Document H21-A2; Vol. 11 No 23.

Registered as
DSM Nutritional Products Ltd
Branch Pentapharm
CH-4002 Basel/Switzerland
P +41 61 706 48 48
F +41 61 319 96 19
pentapharm@pentapharm.com
www.pentapharm.com

Intended Use and Application

Calibrator plasmas to calculate LMWH anticoagulant activity in a blood sample based on clotting times determined with Pefakit® PiCT® (Prothrombinase induced Clotting Time, REF 505-01).

Introduction

Pefakit® PiCT® is a plasma based functional assay for the determination of anticoagulant activity based on FXa and/or FIIa inhibition [1]. Anticoagulant activity of low molecular weight heparin (LMWH) in a blood sample can be read from a standard curve established with this kit. It contains plasmas with specified LMWH anticoagulant activities calibrated against WHO LMWH standard [2].

Use of Pefakit® PiCT® Calibrators LMWH

This kit contains vials of plasma with no LMWH anticoagulant activity (CAL1, 0 IU/ml) and vials of plasma with high LMWH anticoagulant activity (CAL2, refer to certificate) corresponding to measurement points 1 and 4. Points 2 and 3 are prepared by dilution of CAL2 1:5 and 1:1.66 using CAL1 as a diluent. CAL2 is calibrated against WHO LMWH standard.

Reagents

Reagent	Content
CAL1	Calibrator LMWH 1 (human plasma) 3 vials (lyophilisate, to be reconstituted in 1.0 ml of deionized water per vial)
CAL2	Calibrator LMWH 2 (human plasma, supplemented with LMWH) 3 vials (lyophilisate, to be reconstituted in 1.0 ml of deionized water per vial)

Incubate reconstituted calibrators in closed vials for 30' at room temperature and swirl gently before use.

Materials required but not provided

- Deionized water
- Calibrated pipettes (50 - 2000 µl)
- Automated or semi-automated coagulation instruments which employ mechanical or optical detection methods

Note: When using automated or semi-automated coagulation analyzers refer always to manufacturer's operator manual or ask for a detailed adaptation protocol.

Storage and Stability

The kit may be used up to the expiry date given on the label when stored unopened at 2-8°C.

Stability of the plasmas after reconstitution (in original packaging/PP tube):

Calibrators	Stability	
CAL1	-20°C	12 months
	15-25°C	12 hours
CAL2	-20°C	12 months
	15-25°C	12 hours


Calibration plasmas after reconstitution should be frozen for storage. Before use, thaw within 3 minutes at 37°C. Freeze only once. Never store them at 2-8°C.

Expected Values

Clotting times may vary when analyzed on different instruments. Minor variations can occur when different batches of the same reagent are used. Typical results with LMWH Calibrator (100% = 1.44 IU/ml) on two different

instruments are shown in the following tables and diagrams:

LMWH Calibration on ACL 9000™		LMWH Calibration on KC4A™ Micro	
Dilution [IU/ml]	Clotting Time [s]	Dilution [IU/ml]	Clotting Time [s]
0.00	17.4	0.00	21.9
0.43	61.9	0.29	56.0
0.86	98.6	0.86	115.5
1.44	139.0	1.44	162.0

Precautions

The calibrators contain products derived from human blood and have therefore to be handled as potentially infectious material.

Bibliography

1. Calatzis A, Spannagl M, Gempeler-Messina P, Kolde HJ, Schramm W, Haas S. The prothrombinase induced clotting test: A new technique for the monitoring of anticoagulants. *Haemostasis* 2000; 30 (Suppl. 2): 172-174
2. National Committee for Clinical Laboratory Standards. Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays. NCCLS Document H21-A2; Vol. 11 No 23.

Registered as
 DSM Nutritional Products Ltd
 Branch Pentapharm
 CH-4002 Basel/Switzerland
 P +41 61 706 48 48
 F +41 61 319 96 19
 pentapharm@pentapharm.com
 www.pentapharm.com