








Factor XI and XII Deficient Plasma

REF 5184004 Factor XI Deficient Plasma, native 5 x 1 mL

REF 5194008 Factor XII Deficient Plasma, native 5 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola

	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / Производител / výrobce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / oppbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potreba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определенний / роčet stanovení / Definicija
	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Destilana Voda		lot / Charge / lotto / lote / lot / sats / serie / Parti / парти / партида номер / šarže / lot / šarže / Serija
	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaksjonsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakčni pufr / Reakcioni pufer		microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator		catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / καταλογος номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδετικό / Конюгат / Конъюгат / Konjugát / Konjugat		ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk/ έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola		stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / растворите или разредете с / zředit anebo rozpustit v / разбавить или растворить в / nafedite nebo rozpustite v / razrediti ili rastvoriti u		substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffer/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer		washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeoppløsningsskoncentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миещ разтвор / Концентрат промывочного раствора / Koncentrát promývacieho roztoku / Koncetrat solucije za ispiranje
	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico in vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рρησ διαγνωστικής εντόβ σωλήνα / за ин vitro диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnosticky prostředek in vitro / in vitro dijagnostika		

PRODUCT DESCRIPTION

INTEND USE

Factor XI and Faktor XII deficient plasmas are used for the determination of the Factor XI or XII by the one-stage method based on the Activated Partial Thromboplastin Time (aPTT).

COMPOSITION

Factor XI- and Faktor XII Deficient Plasmas are lyophilised stable human deficient plasmas with coagulation activities $\leq 3\%$ of the corresponding factors.

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes	- Solutions/buffers: Distilled water	
REF 5410010	Imidazole buffer	50 mL
REF 5277015	CaCl ₂ 25 mmol/l Solution	100 mL
- Reagents		
REF 5035105	Siron LS	2 x 4 mL
REF 5035118	Siron LIS	2 x 4 mL
REF 5035060	Daptin TC	5 x 2 mL
- Calibrators and Control Plasmas ***		
REF 5220110	Coagulation Reference	5 x 1 mL
REF 5021055	Coagulation Control A	5 x 1 mL
REF 5020040	Coagulation Control N	5 x 1 mL

*** or any other package sizes, special Ceveron® alpha or Technoclot Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- This lot of reagents prepared from human blood and each single plasma used for this lot are HB_sAg, HIV 1/2 Ab and HCV Ab negative.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...+8 °C.

Stability after reconstitution:

RT*	-20 °C
2 hours	1 month

Upon storage, caps should be screwed tightly. * = room temperature

When storing, the tubes must be securely capped.

Deep frozen reagent must be thawed for at least 10 minutes at 37°C and mixed thoroughly before use.

Repeated freezing is not recommended

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of Sodium Citrate solution (0.11 mol/L) and centrifuge for 15 min at a RCF of at least 2500 (corresponding to DIN 58905). Store the plasma at room temperature (up to 4 hours). Before carrying out the test, the patients plasma is diluted 1:10 (0.05 mL + 0.45 mL) using the imidazole buffer. For very low and very high Factor XI / XII levels, however, other dilutions should be used. If the presence of inhibitors is suspected, various dilutions of the plasma sample should be tested.

PREPARATION OF REAGENT

Reconstitute the lyophilised reagents in the prescribed quantities of distilled water and allow them to stand for 10 min at room temperature.¹ Prewarm the CaCl₂ solution to 37°C.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides Application sheets for Ceveron® alpha. The Application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the Application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Pipetting scheme:

0,10 mL	Factor XI or XII Deficient Plasma
+ 0,10 mL	diluted plasma sample
+ 0,10 mL	PTT reagent (shake briefly and incubate for 4 min. at 37 °C)
+ 0,10 mL	CaCl ₂ 25 mmol/l solution 37°C
	determine the point of coagulation

ANALYSES RESULTS

REFERENCE RANGE

Factor XI: 70-120 % of normal, Factor XII: 70-150% of normal

CALCULATION OF THE RESULTS

The predilution 1:10 of the plasma samples is not to be considered in the evaluation.

The plasma samples prediluted at a ratio of 1:10 may be read off directly from the calibration curve. If dilution ratios other than 1:10 are used, the % F XI / XII values of the calibration curve may be converted by using the following formula:

$$\frac{\% \text{ FXI / FXII content}}{10} \times \text{dilution} = \% \text{ FXI / FXII of the sample}$$

CALIBRATION CURVE

Reconstitute the Coagulation Reference. A predilution of 1:10 is attained by using an imidazole buffer in the ratio 1:10 (1 part Coagulation Reference plus 9 parts buffer). Prepare a geometric series of dilutions (1:1 to 1:32) of the predilution (1:10). The 1:1 indication corresponds to the predilution 1:10.

Predilution 1:10	Dilutions of calibration curve					
	1:1	1:2	1:4	1:8	1:16	1:32
% FXI, FXII	100	50	25	12,5	6,3	3,1

Determine the coagulation times of the different dilutions.

Factor XI determination: plot on semi-log paper (x-axis: log activity in %; y-axis: coagulation time in sec)

Factor XII determination; plot on log-log paper (x-axis: log activity in %; y-axis: log coagulation time in sec)

QUALITY CONTROL

In order to verify the accuracy of the results, Coagulation Control A and N should always be tested in the same way as the patient's plasma.

LIMITATION OF THE TEST

Incorrect sample handling can lead to partial activation of the coagulation factors and to falsely elevated single factor determinations.

Lupus anticoagulant can affect the apparent factor activity in single factor determination.

A new calibration is required for each batch of reagents where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced after following a major service of either instruments or equipment.

APPLICATIONS FOR INSTRUMENTS

Application Sheets are available from Technoclone or your local distributor upon request

LITERATURE

Please contact Technoclone or your local distributor.

¹ For standardisation a reconstitution time of 30 min is recommended.

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