|                             |  |                          |
|-----------------------------|--|--------------------------|
| <b>REF 5450701</b>          | <b>TECHNOZYM ADAMTS13 Activity ELISA 96 T.</b>         | <b>ENG</b><br><b>DEU</b> |
| <b>REF 5450761</b>          | <b>TECHNOZYM ADAMTS13 Activity CAL Set 6 x 0.5 mL</b>  |                          |
| <b>REF 5450763</b>          | <b>TECHNOZYM ADAMTS13 Activity CONT Set 2 x 0.5 mL</b> |                          |
| 3050192A Rev.020 16/06/2020 |  | <b>IVD</b> <b>CE</b>     |

## TECHNOZYM ADAMTS13 Activity ELISA - English

### INTENDED USE

The TECHNOZYM ADAMTS13 Activity ELISA is a quantitative test for aid in the diagnosis of TTP, using venous-drawn fresh and/or frozen human citrated (3.2 % sodium citrate) platelet poor plasma, only when considered in conjunction with other clinical indicators. The chromogenic assay is performed on microplate readers capable of reading a wavelength of 450 nm. ADAMTS13 activity measurement serves as an aid in the diagnosis of TTP. The TECHNOZYM ADAMTS13 Activity ELISA is intended for prescription use in laboratories by professionals, qualified to perform ELISA-based assays.

The TECHNOZYM ADAMTS13 Activity CAL Set is used for calibrating the TECHNOZYM ADAMTS13 Activity ELISA.

The TECHNOZYM ADAMTS13 Activity CONT Set is used for quality control of the TECHNOZYM ADAMTS13 Activity ELISA.

### SUMMARY

ADAMTS13 (a disintegrin-like and metalloproteinase with thrombospondin type 1 motif 13) is an enzyme (wWF-cleaving protease or wWF-CP) that specifically cleaves von Willebrand factor (wWF), which induce platelet thrombus formation under high shear stress. If the activity of ADAMTS13 is lowered for some reason, however, unusually large wWF multimers may accumulate, causing thrombosis due to platelet aggregation, which in turn may lead to TMA (thrombotic microangiopathy) such as TTP (thrombotic thrombocytopenic purpura).

### REAGENTS

The TECHNOZYM ADAMTS13 Activity ELISA. contains:

|                | Reagent / Content          | Description  |
|----------------|----------------------------|--|
| 12 x 8 wells   | ELISA test stripes         | Microtiterplate coated with a monoclonal anti GST-Antibody. Drying agent included in aluminum bag    |
| 2 x 6 mL       | GST-wWF73 Substrate        | lyophilized  |
| 6 x 1 x 0.5 mL | Calibrator plasma*         | Numbered from 1 to 6; lyophilized, with lot-specific concentrations (values see batch table)         |
| 2 x 1 x 0.5 mL | Control plasma**           | High and low control plasma, lyophilized, with lot- specific concentrations (values see batch table) |
| 1 x 30 mL      | Reaction buffer            | Liquid, ready to use   |
| 1 x 12 mL      | Conjugate                  | HRP conjugated monoclonal Antibody; liquid, ready to use   |
| 1 x 12 mL      | TMB Substrate              | Tetramethylbenzidine substrate; liquid, ready to use   |
| 1 x 53 mL      | Wash buffer concentrate    | 10-fold concentrated   |
| 1 x 12 mL      | Stop solution              | Sulphuric acid 0.5 M; liquid, ready to use   |
| 1 x 1 plate    | Sample dilution microplate | Only for Sample dilution! (Not coated)   |
| 3 pcs          | Plate sealer               | ELISA Plate sealer during incubation process   |

\*Additionally, the calibrators are sold separately as TECHNOZYM ADAMTS13 Activity CAL Set 6 x 0,5 mL (REF 5450761).

\*\*Additionally, the controls are sold separately as TECHNOZYM ADAMTS13 Activity CONT Set 2 x 0,5 mL (REF 5450763).

#### Material required (not supplied with the kit)

- Distilled water
- Measuring cylinder (500 mL)
- Precision pipettes (5, 50, 100 and 1000 µL)
- Variable pipette (200 and 1000 µL)
- Multichannel and/or dispensing pipettes (100 and 200 µL)
- ELISA washer or multichannel pipette
- ELISA reader with 450 nm filter, with a 620 nm reference filter if available
- Laboratory timer

### Warning and precautions

- IVD for in vitro diagnostic use only.
- This kit is intended for use by personnel trained in laboratory procedures and universal precautions for the use of chemicals and potentially biohazardous substances must be applied.
- All human blood or plasma products as well as test samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas are made from human blood and any individual plasma involved in the procedure is HbsAg, HIV 1/2 Ab and HCV-Ab-negative as tested by FDA approved methods. However, all human blood products should be handled as potentially infectious material.
- Get a Material Safety Data Sheet for this product from www.technoclone.com.

| Symbol   | Warning and Precautions   | Product       |
|--|---|---------------|
|  | H315 causes skin irritation <p>P264 wash hands thoroughly after handling</p> <p>Contains sulphuric acid</p> | Stop solution |

#### Stability and storage

The expiry date printed on the labels is only applicable to storage of the unopened containers at 2...8 °C.

Stability opened / in use:

| Material / Reagent                | State                | Storage   | Stability   |
|-----------------------------------|----------------------|---|-------------|
| ELISA test strips                 | After opening        | 2...8 °C with adhesive film in aluminiumbag with drying agent | Expiry date |
| GST-wWF73 Substrate               | After reconstitution | ≤ -20 °C  | 6 weeks     |
| Calibrators, control plasmas      | After reconstitution | ≤ -20 °C  | 6 months    |
| Reaction Buffer                   | After opening        | 2...8 °C  | 6 months    |
| Conjugate                         | After opening        | 2...8 °C  | 4 months    |
| TMB Substrate                     | After opening        | 2...8 °C  | Expiry date |
| Wash Buffer (10-fold concentrate) | After opening        | 2...8 °C  | 6 months    |

|                |                             |          |             |
|----------------|-----------------------------|----------|-------------|
| Washing Buffer | 1+9 dilution of concentrate | 2...8 °C | 3 weeks     |
| Stop Solution  | after opening               | 2...8 °C | Expiry date |

### TEST PROCEDURE

#### Preparation of plasma samples

Collect nine parts of freshly drawn venous blood in one part trisodium citrate (3.2 %). Refer to CLSI Document H21-A5 for instructions on specimen collection, handling, and storage.

Thaw frozen samples rapidly at 37 °C and centrifuge if necessary. Gently mix before testing. After thawing, the assay must be performed within 2 hours. Samples may be frozen once at -20 °C.

When freezing samples the minimum volume should be 150 µL!

- Sample dilution:* Dilute samples, 31-fold with reaction buffer in a sample dilution Microplate.

Example: 150 µL reaction buffer + 5 µL sample....

For higher precision, volumes can be up scaled, using larger tubes for dilution:

e.g. 600 µL reaction buffer + 20 µL sample

#### Preparation of reagents

Before starting the test, all the required components are to be brought to room temperature.

When reconstituting plasmas, mixing reagents or buffers avoid foaming.

- Washing buffer:* Dilute 1 part by volume washing buffer concentrate with 9 parts by volume distilled water (1+9). Mix well! (diluted washing buffer concentrate = washing buffer). There may be crystalline precipitations which will dissolve at 37 °C within 10 minutes.
- GST-wWF73 Substrate solution:* Substrate solution is reconstituted with 6 mL distilled water and mixed for 10 seconds after a reconstitution time of 15 minutes.
- Calibrators and control plasmas:* Calibrators and control plasmas are reconstituted with 500 µL distilled water and mixed for 10 seconds, after a reconstitution time of 15 minutes. Dilute calibrators and controls 31-fold with reaction buffer in a sample dilution Microplate. (See sample dilution) Reconstituted components are clear to slightly turbid. When freezing calibrator or control plasmas the minimum volume should be 150 µL!

#### Performance of the test

|  |   |                                      |
|--|---|--------------------------------------|
| GST-wWF73 SUBSTRATE INCUBATION (reference 1,2,3,7,9) | Add GST-wWF73 Substrate solution to anti-GST coated test strips <p>Incubate at room temperature</p>       | 100 µL <p>60 minutes</p>             |
| WASHING (reference 1,3,4)                            | Washing buffer  | 3 x 300 µL                           |
| SAMPLE INCUBATION (reference 1,2,5,6,7,9,10)         | Pipette diluted calibrators, control plasmas, samples into test wells <p>Incubate at room temperature</p> | 100 µL <p>30 minutes</p>             |
| WASHING (reference 1,3,4)                            | Washing buffer  | 3 x 300 µL                           |
| CONJUGATE REACTION (reference 1,2,7,9)               | Pipette conjugate into wells <p>Incubate at room temperature</p>  | 100 µL <p>60 minutes</p>             |
| WASHING (reference 1,3,4)                            | Washing buffer  | 3 x 300 µL                           |
| TMB COLOR REAGENT REACTION (reference 1,2,7,9)       | Pipette TMB substrate into test wells <p>Incubate at room temperature</p>                                 | 100 µL <p>30 minutes</p>             |
| STOPPING (reference 1,2,7)                           | Pipette stopping solution into wells  | 100 µL                               |
| MEASUREMENT (reference 8)                            | ELISA reader, 450 nm  | shake 5 sec., measure within 10 min. |

### References

- Reagents of different lots must not be combined.
- Precision and performance, among others, essentially depend on the following factors:
  - Thorough mixing of all substances used for dilution, 10 sec. with vortex mixer.
  - Run calibrators, controls and samples in duplicates.
  - Incubate at indicated temperature (room temperature 18...25 °C).
  - Strict observance of the order of pipetting and of the time element as indicated.
  - The time for sample incubation, conjugate and substrate reaction as indicated, starts after pipetting the last sample. Incubation times should not vary by more than ± 5 %.
  - During sample incubation and conjugate reaction, the time for pipetting calibrators / control plasmas / samples and / or conjugate solutions must not exceed 60 seconds per ELISA test strip (8 wells).
  - During substrate reactions and at stopping, the time needed for pipetting the substrate and / or the stopping solution must not exceed 10 seconds per ELISA test strip. Short pipetting times may be secured by using multichannel pipettes.
- Label / number strips with a water-resistant pen, in case the strips accidentally fall out of the frame during testing
- After the last washing, wells must be aspirated thoroughly, turned upside down and positioned on a blotting paper. The last remnants must be removed by gentle tapping.
- A calibration curve has to be created for every assay.
- Samples / calibrators / controls can be transferred from sample dilution microplate to anti-GST microplate by multichannel pipette. Do not forget to change tips for every strip!
- No agitation is required during each reaction step.
- By measuring the difference in wavelength at 450 nm and 620 nm, the precision of the test is increased.
- For every incubation step, the test plate has to be covered with plate sealer.
- It is not mandatory to use all six calibrators for creating a calibration curve. For sample screening, a calibration curve using calibrators 1, 2, 3, 4, and 6 is sufficient. When the focus is particularly in the lower range of ADAMTS13 activity, a calibration curve with calibrators 2, 3, 4, 5 and 6 is adequate.

### LIMITATION OF THE TEST

Samples containing EDTA cannot be used, because EDTA is a strong inhibitor of ADAMTS13 function.

Hemolysis: No interference is observed with samples containing up to 200 mg/dL hemoglobin, which corresponds with a moderate hemolysis.

Lipemia: No interference is observed with samples containing up to 300 mg/dL Intralipid™, which corresponds with a moderate to severe concentration.

Icterus: No interference was observed with samples containing up to 15 mg/dL bilirubin (conjugated as well as unconjugated), which corresponds with a moderate to severe concentration.

Rheumatoid factor: No interference was observed up to 28 IU/mL, with corresponds with a 2-fold concentration of normal.

Anti CD20 antibodies: No interference was observed up to a level of 200 µg/mL, which corresponds to the upper level of serum concentrations found after Rituximab administration.

### INTERPRETATION OF RESULTS

TECHNOZYM ADAMTS13 Activity results are reported in IU/mL. The assay results should be used with other information, including the clinical context, in forming a diagnosis.

### CALCULATION OF RESULTS

#### Setting up a reference curve:

X axis: ADAMTS13 Activity [IU/mL]

Y axis: Extinction at 450 nm

Graph plot is linear-linear with a best fit.

#### Assessment of reference curve:

The validity of the test may be checked based on the calculated control values.

#### Measuring concentration of samples:

Read off the concentration from the reference curve.

If there are samples, with extinction coefficients higher than the extinction of the highest point on the calibration curve, they have to be pre diluted with reaction buffer (1+1 or 1+3). The measured concentration then has to be multiplied with the dilution factor 2 or 4, respectively.

### REFERENCE RANGE

Normal range for ADAMTS13 Activity: 0.4 – 1.3 IU/mL.

It is recommended that individual laboratories establish their own normal range. When interpreting the serological results, the history of the patient has to be taken into account.

### PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

### Clinical Performance

Performance was determined in patients with proven TTP.

| Parameter                 | Value                  | 95 <span> </span> % CI                           |
|---------------------------|------------------------|--|
| Sensitivity               | 98.73 <span> </span> % | 93.15 <span> </span> % to 99.97 <span> </span> % |
| Specificity               | 92.86 <span> </span> % | 80.52 <span> </span> % to 98.50 <span> </span> % |
| Positive Predictive Value | 96.30 <span> </span> % | 89.73 <span> </span> % to 98.72 <span> </span> % |
| Negative Predictive Value | 97.50 <span> </span> % | 84.74 <span> </span> % to 99.64 <span> </span> % |

Patients with suspected and clinically proven Thrombotic Thrombocytopenic Purpura (TTP) were tested. Based on a 0.1 IU/mL cut off to diagnose TTP, assay performance parameters were calculated.

### Precision

Reproducibility was determined with different samples.

| Sample code  | Assigned value [IU/mL] | CV <span> </span> % within run | CV <span> </span> % total |
|--------------|------------------------|--------------------------------|---------------------------|
| High control | 0.909                  | 7.7                            | 10.5                      |
| Low control  | 0.239                  | 9.3                            | 13.8                      |

#### Limit of quantification and assay range

When assay is performed as indicated in section 'Test Procedure', the limit of quantification of this assay is 0,0071 IU/mL ADAMTS13 activity. The upper limit of the assay range may vary with each lot of kit depending on the assayed value of the calibrator plasma supplied in the kit. Samples with values outside the range of the reference curve should be re-tested at an appropriate dilution to obtain accurate results.

### STANDARDISATION

Standards and controls are produced from a normal donor and they are calibrated against the WHO International Standard for ADAMTS13. Consult the batch table.

### LITERATURE

Please contact Technoclone www.technoclone.com or your local distributor.

### EDITORIAL NOTE

This document is available in several languages. The translations have been done using the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

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**TECHNOZYM** is a registered trademark of Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH.
TECHNOZYM is a registered trademark of Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH.
Technoclone Herstellung von Dignostika und Arzneimitteln GmbH, Brunner Str. 67 - 1230 Vienna, Austria