NOTICE TECHNIQUE

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XI. BIBLIOGRAPHIE

1. Biggs R. Human blood coagulation, haemostasis and thrombosis 3rd ed. Oxford: Blackwell Scientific Publications; 1984.
2. Henry, John Bernard, MD, (ed) Clinical Diagnosis and Management by laboratory Methods, Philadelphia, PA, W. B. Saunders Co., 1984, pp. 51-92
3. Triplett DA, Smith C. Routine testing in the coagulation laboratory. In: Triplett DA, editor. Laboratory evaluation of coagulation. Illinois: ASCP Press; 1982. p. 28-51.
4. Biosafety in microbiological and biomedical laboratories, 4th ed. Centers for Disease Control and Prevention/ National Institutes of Health, 1999
5. Collection, transport, and processing of blood specimens for coagulation testing and general performance of coagulation assays, NCCLS, H21-A3. 1998.
6. Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989.
7. CLIA 1988 – Code of Federal Regulations, 42CFR493.1253, 1998.
8. Determination of Factor Coagulant Activities, NCCLS, H48-A. 1997.
9. Gilmer PR. Preanalytical variables in coagulation testing. In: Triplett DA, editor. Laboratory evaluation of coagulation. Illinois: ASCP Press; 1982. p. 1-8.

XII. FABRICANT

CRYOPEP
ZAC Parc 2000
83 rue Yves Montand
34080 Montpellier
Tel : 04 67 10 71 20
Fax : 04 67 10 71 21
contact@cryopep.com
www.cryopep.fr

6-1800-05 Factor VIII Inhibitor Plasma - Weak Control

I. INTENDED USE

Plasma must be used according to the procedures described in the instructions of the reagents.

For In Vitro Diagnostic

II. SUMMARY AND PRINCIPLE

Deficiencies in coagulation factors may have congenital or acquired etiologies and can compromise in vivo hemostasis¹.

These plasmas can be used for clinical laboratory use in conjunction with a Factor VIII activity assay as weak positive control during the determination of the factor VIII inhibitors using Bethesda assay or a modified Nijmegen-Bethesda assay^{2,3}.

The inhibitory capacity is indicated in Bethesda Units / mL on the certificate of analysis.

III. PRESENTATION

Kit - 25 vials of 0.5 mL.

IV. REAGENTS

Factor VIII Inhibitor Plasma is produced from native Factor VIII Deficient human plasma to which an inhibitor antibody was added in order to provide a precise neutralizing activity.

ATTENTION : All blood products should be treated as potentially infectious. Source material from which this product was derived was found to be negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. Accordingly, these human blood based products should be handled and discarded as recommended for any potentially infectious human specimen. According to the national institute of health regulations and the recommendations of level 2 from the manual of biosafety in microbiological and biomedical laboratories⁴, 1999.

Storage

Once received in dry ice, place the plasma in the freezer so that it stabilizes. When stored at -40 to -80°C, this plasma is stable to the end of the month indicated on the product packaging. The product stability is 7 days at -20°C.

Thawed material should not be refrozen

Handling

Thaw each vial at 37°C (± 1°C) in a water bath during 3 min. Thawing times are important and should be strictly respected. Allow thawed plasma to reach room temperature (18 to 25°C) during 15 to 30 min. and invert gently prior to use.

Plasma may be used for up to 2 hours after thawing, if capped in the original vial and maintained at room temperature or up to 4 hours if maintained at 2 to 8°C.

V. INSTRUMENTS

Each laboratory should prepare the local instrument in accordance with the manufacturer's instructions for use.

VI. PROCEDURE

After thawing and preparing deficient plasma, use in accordance with established laboratory procedures for APTT based determination of inhibitors of coagulation factors.

Material provided

Factor VIII native deficient plasma containing inhibitors of Factor VIII.

Materials required but not provided

- Water bath (37 ± 1°C)
- Assay reagents
- Owren's Veronal Buffer or equivalent
- Coagulation instrument or assay system
- Calibration plasma
- Quality control material
- Plastic test tubes
- Timer

Specimen Collection and Preparation

Patient samples should be collected into 105-109 mmol/L sodium citrate dihydrate anticoagulant (3.2%) in a ratio of 9 parts blood to 1 part anticoagulant. Patient plasma is derived by centrifugation at 1500 g for 15 min. and should be tested within four hours of collection when maintained at 2 to 4°C in accordance with NCCLS guidelines⁵.

Quality Control

Each laboratory should establish its own quality control (QC) ranges using acceptable statistical methods. These QC ranges may then be used to monitor and validate the integrity of the test system⁶. For all coagulation tests, the laboratory must include at least two levels of control for every eight hours of operation and any time a change in reagents occurs⁷.

VII. RESULTS

Factor activity values recovered below the normal range may be indicative of a factor deficiency (congenital or acquired). Each laboratory should establish its own normal range for Factor VIII activity in accordance with NCCLS guidelines⁸.

VIII. LIMITATIONS OF THE PROCEDURE

When proper control values are not obtained, assessment of each component of the test system including reagents, control plasmas, instrumentation and operator technique must be undertaken in order to ascertain that all other components are functioning properly⁹.

IX. EXPECTED VALUES

Expected values may vary according to reagent, instrument and technique employed. It is recommended each laboratory establish its own normal range for factor activity.

X. PERFORMANCE CHARACTERISTICS

When used according to recommended methods, results are subject to the limitations of the assay system (i.e. reagents, instrument) in use.

XI. BIBLIOGRAPHY

1. Biggs R. *Human blood coagulation, haemostasis and thrombosis* 3rd ed. Oxford: Blackwell Scientific Publications, 1984.
2. Henry John Bernard (MD), (ed) *Clinical Diagnosis and Management by Laboratory Methods*, Philadelphia, PA, W. B. Saunders Co., 1984, pp. 51-92
3. Triplett DA, Smith C. Routine testing in the coagulation laboratory. In: Triplett DA, editor. *Laboratory evaluation of coagulation*. Illinois: ASCP Press; 1982. p. 28-51.
4. *Biosafety in microbiological and biomedical laboratories*, 4th ed. Centers for Disease Control and Prevention/ National Institutes of Health, 1999
5. Collection, transport, and processing of blood specimens for coagulation testing and general performance of coagulation assays, NCCLS, H21-A3. 1998.
6. Cembrowski GS, Carey RN. *Laboratory quality management*. Chicago: ASCP Press; 1989.
7. CLIA 1988 – *Code of Federal Regulations*, 42CFR493.1253, 1998.
8. Determination of Factor Coagulant Activities, NCCLS, H48-A. 1997.
9. Gilmer PR. Preanalytical variables in coagulation testing. In: Triplett DA, editor. *Laboratory evaluation of coagulation*. Illinois: ASCP Press; 1982. p. 1-8.

XII. MANUFACTURER

CRYOPEP
ZAC Parc 2000
83 rue Yves Montand
34080 Montpellier
Tel : 04 67 10 71 20
Fax : 04 67 10 71 21
contact@cryopep.com
www.cryopep.fr