
NOTICE TECHNIQUE

PLASMAS DEFICIENTS NATIFS

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X. BIBLIOGRAPHIE

1. Biggs R. Human blood coagulation, haemostasis and thrombosis 3rd ed. Oxford: Blackwell Scientific Publications; 1984.
2. Henry, John Bernard, MD, (ed) Clinical Diagnosis and Management by laboratory Methods, Philadelphia, PA, W. B. Saunders Co., 1984, pp. 51-92
3. Triplett DA, Smith C. Routine testing in the coagulation laboratory. In: Triplett DA, editor. Laboratory evaluation of coagulation. Illinois: ASCP Press; 1982. p. 28-51.
4. Biosafety in microbiological and biomedical laboratories, 4th ed. Centers for Disease Control and Prevention/ National Institutes of Health, 1999
5. Collection, transport, and processing of blood specimens for coagulation testing and general performance of coagulation assays, CLSI, H21-A3. 1998.
6. Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989.
7. CLIA 1988 – Code of Federal Regulations, 42CFR493.1253, 1998.
8. Determination of Factor Coagulant Activities, CLSI, H48-A. 1997.
9. Gilmer PR. Preanalytical variables in coagulation testing. In: Triplett DA, editor. Laboratory evaluation of coagulation. Illinois: ASCP Press; 1982:1-8.

XI. FABRICANT

CRYOPEP
ZAC Parc 2000
83 rue Yves Montand
34080 Montpellier
Tel : 04 67 10 71 20
Fax : 04 67 10 71 21
contact@cryopep.com
www.cryopep.fr

TECHNICAL MANUAL

NATIVE DEFICIENT PLASMAS

RUO

-80°C / -60°C

Reference	Native Deficient Plasma
7-0200	Factor II Deficient
7-0500	Factor V Deficient
7-0700	Factor VII Deficient
7-0800	Factor VIII Deficient
7-0900	Factor IX Deficient
7-1000	Factor X Deficient
7-1100	Factor XI Deficient
7-1200	Factor XII Deficient
7-1300	Factor XIII Deficient
7-1700	Prekallikrein Deficient
7-1800	Factor VIII Deficient with inhibitor
7-2000	Protein C Deficient

I. INTENDED USE

These native deficient plasmas are recommended to evaluate the activity of clotting factors by the determination of prothrombin time (PT) or partial thromboplastin time activator (APTT) requiring the use of a plasma deficient in factor (< 1 %). Plasmas should be used in these procedures according to reagent insert directions.

II. SUMMARY AND PRINCIPLE

The clotting factor deficiencies may have congenital or acquired origin and can undermine the process of in vivo hemostasis¹. When doing coagulation assays, Good Laboratory Practices require that procedures be monitored with a suitable control product^{2,3}. These products may be used in coagulation assays using citrated human plasmas.

III. REAGENTS

Coagulation factor deficient plasmas: Fresh frozen citrated human plasma from congenital coagulation factor deficient donors.

WARNINGS: Each individual donor plasma used in these products has been tested by an FDA approved method and found to be non-reactive for the presence of HBsAg and antibody to HIV. Because no known test method can offer complete assurance that these or other infectious agents are absent, this product should be handled at the Biosafety Level 2 as recommended for any human blood-based product in the Centers for Disease Control / National Institutes of Health manual "Biosafety in Microbiological and Biomedical laboratories", 1999⁴.

Storage

Recommended storage: -40°C/-80°C until to expiry date.
Plasmas must remain frozen.
May be stored at -20°C if used within 7 days.

DO NOT REFREEZE

Preparation

Place vial in 37±1°C water bath until plasma is thawed. Exact time is determined by volume of plasma in vial.

TABLE OF DEFROST	
Volume of plasma	Water-bath at 37±1°C
1.0 mL	4 min
0.5 mL	3 min

Mix gently and leave plasma stabilize at room temperature (18 à 25°C). Plasma should be used within 2 hours after being thawed at room temperature or within 4 hours if stored at +2/+8°C. Plasma must be discarded once thawed and used.

IV. INSTRUMENTS

Each laboratory must prepare the necessary instruments stipulated in the manufacturer's instructions.

V. PROCEDURE

After defrosting and preparing the deficient plasma, use the plasma, as described, according to the procedures established in the laboratory for quantitative chronometric assays of the coagulation factors.

Provided material

Native deficient plasma

Materials requested but not provided

- Water bath (37 ± 1 °C)
- Reagents dosing
- Buffer Owren -Koller or equivalent
- Instrument or coagulation assay system
- Calibration plasma
- Material quality control
- Tubes and plastic pipettes
- Stopwatch
- Paper Log-Log

Sample collection & Preparation

The blood samples of patients should be collected on anticoagulant (trisodium citrate solution 3.2 % - 109 mmol / L) in a ratio of 9 parts of blood to 1 part of anticoagulant. The patient's plasma is obtained by centrifugation at 1500 g for 15 minutes and should be tested in 4 hours after collection when it is kept at +2/+4°C as stated in the instructions CLSI⁵.

Quality Control

Each laboratory should establish its own quality control standards using acceptable statistical methods. These standards should be used to control and validate the integrity of test⁶. For all coagulation tests, the laboratory must include at least two levels of control every 8 hours and in any case, a change of reagents should be done⁷.

VI. RESULTS

Coagulation factor activity values found below normal may indicate factor deficiency (congenital or acquired). Each laboratory must establish its own range for the activity of coagulation factors in accordance with the instructions of the CLSI⁸.

VII. LIMITES DE LA METHODE

When the expected values of the controls do not comply, the control of each component of the measurement system (reagents, control plasmas, operating instruments and methods) must be carried out to ensure that all the components are functionally correct⁹.

VIII. EXPECTED VALUES

Expected values are dependent on instrumentation, reagent system and methodology used. Due to the wide variation in these components, it is recommended that these plasmas be evaluated and quality controlled by the user for their specific system.

IX. PERFORMANCE CHARACTERISTICS

As with Expected Values, the performance characteristics are dependent on instrumentation, reagents and specific methodologies employed to do the testing.

X. BIBLIOGRAPHY

1. Biggs R. Human blood coagulation, haemostasis and thrombosis 3rd ed. Oxford: Blackwell Scientific Publications; 1984.
2. Henry, John Bernard, MD, (ed) Clinical Diagnosis and Management by laboratory Methods, Philadelphia, PA, W. B. Saunders Co., 1984, pp. 51-92
3. Triplett DA, Smith C. Routine testing in the coagulation laboratory. In: Triplett DA, editor. Laboratory evaluation of coagulation. Illinois: ASCP Press; 1982. p. 28-51.
4. Biosafety in microbiological and biomedical laboratories, 4th ed. Centers for Disease Control and Prevention/ National Institutes of Health, 1999
5. Collection, transport, and processing of blood specimens for coagulation testing and general performance of coagulation assays, CLSI, H21-A3. 1998.
6. Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989.
7. CLIA 1988 – Code of Federal Regulations, 42CFR493.1253, 1998.
8. Determination of Factor Coagulant Activities, CLSI, H48-A. 1997.
9. Gilmer PR. Preanalytical variables in coagulation testing. In: Triplett DA, editor. Laboratory evaluation of coagulation. Illinois: ASCP Press; 1982:1-8.

XI. MANUFACTURER

CRYOPEP
ZAC Parc 2000
83 rue Yves Montand
34080 Montpellier
Tel : 04 67 10 71 20
Fax : 04 67 10 71 21
contact@cryopep.com
www.cryopep.fr